

*Evacuated blood collection
made easy*

NEW



Helping all people
live healthy lives

BD Vacutainer® Flashback Needle

Venous blood collection
with greater confidence

A Unique Design

For greater comfort and confidence in use

BD Vacutainer® Flashback needle gives you instantaneous confirmation of the needle placement, adding greater comfort in use.

- Avoids needle probing and the multiplication of attempts to draw blood from a patient.
- Ensures easy training of new staff or inexperienced personnel, and full utilization of closed-system collection.



Now, the technique of closed-system blood collection is easier for everyone.

Reliable results start with reliable samples

Closed blood collection system:

- Better specimen quality, less needs for redraws
- Less manipulations
- Less exposure to blood
- Faster procedures

BD Vacutainer® Flashback needle

Ref.	Gauge	Length	Colour Coding	Packaging
301746	21G – 0.8 mm	1" – 25 mm		50 /box 1000/case
301747	22G – 0.6 mm	1" – 25 mm		50 /box 1000/case
301748	22G – 0.6 mm	1,25" – 32 mm		50 /box 1000/case

BD Vacutainer® Flashback needle ideally complements our range of blood collection devices and accessories. Ask your BD representative for more information.



BD Diagnostics
Preanalytical Systems
21 Between Towns Road, Cowley, OX4 3LY
Tel: 01865 748844 Fax: 01865 781528
www.bd.com